

Please amend the above-referenced application as follows:

**In The Specification:**

Please replace the title at page 1, lines 1 and 2, with the following re-written title:

B. Apolipoprotein Biopolymer Markers [Predictive] Indicative Of  
Type II Diabetes

Please replace the paragraph beginning at page 40, line 12, with the following rewritten paragraph:

Preparatory Protocols:

Any of these protocols may be selected from a column flow-through stream, a column elution stream, or a column scrub stream.

Hi Q is a strong anion exchanger made of methyl acrylate co-polymer with the functional group:  $-N^+(CH_3)_2$ ;

Hi S is a strong cation exchanger made of methyl acrylate co-polymer with the functional group:  $-SO_3^-$ ;

b<sub>1</sub> DEAE is a diethylaminoethyl which is a weak cation exchanger made of methyl acrylate co-polymer with the functional group  $-N^+(C_2H_5)_2$ ;

PS is phenyl [sepharose] SEPHAROSE;

BS is buytl [sepharose] SEPHAROSE.

Please replace the paragraph beginning at page 41, line 2,  
with the following rewritten paragraph:

---

b2 Note that the supports, i.e. methyl acrylate and [sepharose]  
SEPHAROSE are different, but non-limiting examples, as the same  
functional group on different supports will function, albeit  
possibly with different effects.

---

Please replace the paragraph beginning at page 41, line 20,  
with the following rewritten paragraph:

---

Butyl [sepharose] SEPHAROSE column protocol:

- b3
- 1) Cast 150  $\mu$ l bed volume column;
  - 2) Equilibrate column in 5 bed volumes of 1.7 M  
( $\text{NH}_4$ )<sub>2</sub>SO<sub>4</sub> in 50 mM PB pH 7.0 (binding buffer);
  - 3) Dissolve 35  $\mu$ l of sera in 465  $\mu$ l of binding buffer  
and apply;
  - 4) Wash column in 5 bed volumes of binding buffer;
  - 5) Elute column in 120  $\mu$ l of 0.4 M ( $\text{NH}_4$ )<sub>2</sub>SO<sub>4</sub> in 50 mM PB  
pH 7.0;
  - 6) Elute column in 120  $\mu$ l of 50 mM PB pH 7.0;
  - 7) Scrub column with 120  $\mu$ l sequentially with each of  
0.1% triton, 1.0% triton and 2% SDS in 62.5 mM Tris pH 6.8.
-

Please replace the paragraph beginning at page 42, line 12,  
with the following rewritten paragraph:

---

Phenyl [sepharose] SEPHAROSE column protocol:

- 1) Cast 150  $\mu$ l bed volume column;
  - 2) Equilibrate column in 5 bed volumes of 1.7 M  $(\text{NH}_4)_2\text{SO}_4$  in 50 mM PB pH 7.0 (binding buffer);
  - 3) Dissolve 35  $\mu$ l of sera in 465  $\mu$ l of binding buffer and apply;
  - 4) Wash column in 5 bed volumes of binding buffer;
  - 5) Elute column in 120  $\mu$ l of 0.2 M  $(\text{NH}_4)_2\text{SO}_4$  in 50 mM PB pH 7.0;
  - 6) Elute column in 120  $\mu$ l of 50 mM PB pH 7.0;
  - 7) Scrub column with 120  $\mu$ l sequentially with each of 0.1% triton, 1.0% triton and 2% SDS in 62.5 mM Tris pH 6.8.
-

Please replace the paragraph beginning at page 66, line 2, with the following rewritten paragraph:

---

*be* The instant invention involves the use of a combination of preparatory steps in conjunction with mass spectroscopy and time-of-flight detection procedures to maximize the diversity of biopolymers which are verifiable within a particular sample. The cohort of biopolymers verified within such a sample is then viewed with reference to their ability to evidence at least one particular disease state; thereby enabling a diagnostician to gain the ability to characterize either the presence or absence of [said] at least one disease state relative to recognition of the presence and/or the absence of [said] the biopolymer, predict disease risk assessment, and develop therapeutic avenues against [said] the disease.

---